

Nur für Forschungszwecke

# MED18 Polyklonaler Antikörper

Katalog-Nr.: 15352-1-AP

1 Publikationen



## Allgemeine Informationen

Katalog-Nr.:

15352-1-AP

Größe:

150ul, Konzentration: 550 µg/ml von Nanodrop und 267 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG7598

GenBank-Zugangsnummer:

BC002694

GeneID (NCBI):

54797

Vollständiger Name:

mediator complex subunit 18

Berechnete Masse:

24 kDa

Beobachtete Masse:

24 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:6000

## Anwendungen

Geprüfte Anwendungen:

WB, ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Maus

Positivkontrollen:

WB: HeLa-Zellen, Jurkat-Zellen

## Hintergrundinformationen

MED18, also named as Mediator of RNA polymerase II transcription subunit 18, is a 208 amino acid protein, which belongs to the Mediator complex subunit 18 family. MED18 is a component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. MED18 as a mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. MED18 is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors.

## Bemerkenswerte Veröffentlichungen

Verfasser

Wen-Shuai Tang

Pubmed ID

34233190

Journal

Cell Rep

Anwendung

WB

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 15352-1-AP (MED18 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.