

Nur für Forschungszwecke

YPEL3 Polyklonaler Antikörper

Katalog-Nr.: 15403-1-AP

Vorgestelltes Produkt

4 Publikationen



Allgemeine Informationen

Katalog-Nr.:
15403-1-AP

Größe:
150ul, Konzentration: 300 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG7626

GenBank-Zugangsnummer:
BC005009

GeneID (NCBI):
83719

Vollständiger Name:
yippee-like 3 (Drosophila)

Berechnete Masse:
14 kDa

Beobachtete Masse:
14 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
IF, IHC, WB

Getestete Reaktivität:
Human, Maus, Ratte

Zitierte Arten:
Human

Hintergrundinformationen

YPEL3 is part of a five-member family of closely related paralogues: YPEL1-5, which are named in reference to their Drosophila orthologue. Studies showed that YPEL3 is regulated by p53 and its encoding protein induces cellular senescence in human tumor and normal cells, indicating that YPEL3 is a p53-dependent, tumor suppressor gene. (PMID:20388804) 15403-1-AP was raised against full length of YPEL3 protein (1-119aa); it can recognize YPEL1-4.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Hyoung-Seok Baek	34675221	Sci Rep	WB, IF
Qing Liu	36911393	Front Genet	IHC
Tuttle Rebecca R	21671470	Int J Cancer	WB

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

