

Nur für Forschungszwecke

SF4 Polyklonaler Antikörper

Katalog-Nr.: 15614-1-AP



Allgemeine Informationen

Katalog-Nr.: 15614-1-AP	GenBank-Zugangsnummer: BC002986	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 650 µg/ml von Nanodrop und 367 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 57794	Empfohlene Verdünnungen: WB 1:500-1:1000 IHC 1:20-1:200 IF 1:10-1:100
Wirt: Kaninchen	Vollständiger Name: splicing factor 4	
Isotyp: IgG	Berechnete Masse: 72 kDa	
Immunogen Katalognummer: AG8006	Beobachtete Masse: 72 kDa	

Anwendungen

Geprüfte Anwendungen: IF, IHC, WB, ELISA	Positivkontrollen: WB : HEK-293-Zellen, Maushodengewebe IHC : humanes Hodengewebe, IF : HepG2-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	
Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HEK-293 cells were subjected to SDS PAGE followed by western blot with 15614-1-AP (SF4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human testis using 15614-1-AP (SF4 antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 15614-1-AP (SF4 antibody) at dilution of 1:100 (under 40x lens).



Immunofluorescent analysis of () fixed HepG2 cells using 15614-1-AP (SF4 antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG.