

Nur für Forschungszwecke

GRB7 Polyklonaler Antikörper

Katalog-Nr.:15738-1-AP



Allgemeine Informationen

Katalog-Nr.: 15738-1-AP	GenBank-Zugangsnummer: BC006535	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 167 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 2886	
Wirt: Kaninchen	Vollständiger Name: growth factor receptor-bound protein 7	
Isotyp: IgG	Berechnete Masse: 532 aa, 60 kDa	
Immunogen Katalognummer: AG8537		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Grb7 belongs to a family of adapter proteins. This family, including Grb10, Grb7 and Grb 14, contains a proline-rich, a RA, a PH, a BPS and SH2 domain. Grb7 associate with a number of receptors including EGFR and FAK. Functions of Grb7 include mobility because it associates receptors such FAK, and direct evidence from overexpression of this proteins. This antibody recognize the C-terminal of GRB7.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

