

Nur für Forschungszwecke

C1QL1 Polyklonaler Antikörper

Katalog-Nr.:15757-1-AP



Allgemeine Informationen

Katalog-Nr.:
15757-1-AP

Größe:
150ul , Konzentration: 400 µg/ml von
Nanodrop und 227 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG8348

GenBank-Zugangsnummer:
BC008798

GeneID (NCBI):
10882

Vollständiger Name:
complement component 1, q
subcomponent-like 1

Berechnete Masse:
258 aa, 26 kDa

Beobachtete Masse:
30 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : Maushirngewebe, Rattenhirngewebe

Hintergrundinformationen

C1QL1 (complement component 1, q subcomponent-like 1) is one of the C1q/tumor-necrosis-factor-related protein (CTRP) family members. C1QL1 is synthesized and secreted mainly by the brain and adipose tissues. It is implicated in synapse formation or elimination in cultured hippocampal neurons in vitro.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 15757-1-AP (C1QL1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.