

TSPAN18 Polyklonaler Antikörper

Katalog-Nr.: 15919-1-AP

Allgemeine Informationen

Katalog-Nr.:

15919-1-AP

Größe:

150ul, Konzentration: 173 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG8724

GenBank-Zugangsnummer:

BC019342

GeneID (NCBI):

90139

Vollständiger Name:

tetraspanin 18

Berechnete Masse:

267 aa, 29 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

TSPAN18 is a membrane protein of the tetraspanin superfamily, which are characterized by the presence of four conserved transmembrane regions. Tetraspanins have been involved in diverse processes such as cell activation and proliferation, adhesion and motility, differentiation, and cancer. TSPAN18 gene mutation has been associated with Schizophrenia in Han Chinese (PMID: 22037552; 23505562).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

