

Nur für Forschungszwecke

FRK Polyklonaler Antikörper

Katalog-Nr.:15936-1-AP



Allgemeine Informationen

Katalog-Nr.:
15936-1-AP

Größe:
150ul, Konzentration: 340 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG10024

GenBank-Zugangsnummer:
BC012916

GeneID (NCBI):
2444

Vollständiger Name:
fyn-related kinase

Berechnete Masse:
505 aa, 56 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

The fyn-related kinase (FRK) belongs to the tyrosine kinase family of protein kinases and is also named as PTK5 (protein-tyrosine kinase 5), RAK (Nuclear tyrosine protein kinase). The Frk is a tyrosine kinase expressed in various tissues, such as the gastrointestinal tract and pancreatic islet cells. It affects pancreatic beta cell number during embryogenesis and promotes beta cell cytotoxic signals in response to streptozotocin (PMID:20157399). It may function as a tumor suppressor (PMID:19345329).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

