

Nur für Forschungszwecke

MSK1 Polyklonaler Antikörper

Katalog-Nr.: 15977-1-AP



Allgemeine Informationen

Katalog-Nr.:

15977-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG8765

GenBank-Zugangsnummer:

BC017187

GeneID (NCBI):

9252

Vollständiger Name:

ribosomal protein S6 kinase, 90kDa,
polypeptide 5

Berechnete Masse:

549aa, 62 kDa; 802aa, 90 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

IF 1:10-1:100

Anwendungen

Geprüfte Anwendungen:

IF, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

IF : MCF-7-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of MCF-7 cells, using RPS6KA5 antibody 15977-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).