

Nur für Forschungszwecke

INTS4 Polyklonaler Antikörper

Katalog-Nr.: 16130-1-AP



Allgemeine Informationen

Katalog-Nr.:

16130-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von
Nanodrop;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG9116

GenBank-Zugangsnummer:

BC009995

GeneID (NCBI):

92105

Vollständiger Name:

integrator complex subunit 4

Berechnete Masse:

963aa, 108 kDa; 505aa, 57 kDa

Beobachtete Masse:

108 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: HeLa-Zellen, HepG2-Zellen, Jurkat-Zellen, MCF-7-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 16130-1-AP (INTS4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.