

Allgemeine Informationen

Katalog-Nr.:

16223-1-AP

Größe:

150ul, Konzentration: 600 µg/ml von Nanodrop und 400 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG9276

GenBank-Zugangsnummer:

BC021666

GeneID (NCBI):

202374

Vollständiger Name:

serine/threonine kinase 32A

Berechnete Masse:

396aa, 46 kDa; 166aa, 20 kDa

Beobachtete Masse:

46 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: HeLa-Zellen, MCF-7-Zellen

Hintergrundinformationen

STK32A (Serine/threonine-protein kinase 32A) is also named as YANK1 and belongs to Ser/Thr protein kinase family. It can affect cell metabolism by protein serine/threonine kinase phosphorylation. The gene encoding STK32A maps to human chromosome 5 and this gene is associated with Cockayne syndrome through the ERCC8 gene and familial adenomatous polyposis through the adenomatous polyposis coli (APC) tumor suppressor gene.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 16223-1-AP (STK32A antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.