

Nur für Forschungszwecke

FOXS1 Polyklonaler Antikörper

Katalog-Nr.:16234-1-AP



Allgemeine Informationen

Katalog-Nr.:
16234-1-AP

Größe:
150ul, Konzentration: 900 µg/ml von
Nanodrop und 453 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG8957

GenBank-Zugangsnummer:
BC013408

GeneID (NCBI):
2307

Vollständiger Name:
forkhead box S1

Berechnete Masse:
330 aa, 35 kDa

Beobachtete Masse:
35 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

WB : Mauslebergewebe, Rattenlebergewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse liver tissue were subjected to SDS PAGE followed by western blot with 16234-1-AP (FOXS1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.