

TUSC4 Polyklonaler Antikörper

Katalog-Nr.: 16319-1-AP

Allgemeine Informationen

Katalog-Nr.:

16319-1-AP

Größe:

150ul, Konzentration: 93 µg/ml durch 10641

die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG9372

GenBank-Zugangsnummer:

BC021984

GeneID (NCBI):

10641

Vollständiger Name:

tumor suppressor candidate 4

Berechnete Masse:

380 aa, 44 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

NPRL2 (NPR2 Like, GATOR1 Complex Subunit) is a Protein Coding gene associated with Epilepsy, Familial Focal. As a component of the GATOR1 complex functions as an inhibitor of the amino acid-sensing branch of the TORC1 pathway. Homologous to yeast nitrogen permease (candidate tumor suppressor) (NM_006545, BC021984) gene encodes a protein of molecular mass of 50 kDa. Data by using overlapping and nested homozygous deletions methods indicate that this protein might be a tumor suppressor involving in the development of lung cancers.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

