

Nur für Forschungszwecke

CTSZ Polyklonaler Antikörper

Katalog-Nr.:16578-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:
16578-1-AP

Größe:
150ul , Konzentration: 400 µg/ml von
Nanodrop und 267 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG9865

GenBank-Zugangsnummer:
BC042168

GeneID (NCBI):
1522

Vollständiger Name:
cathepsin Z

Berechnete Masse:
303 aa, 34 kDa

Beobachtete Masse:
34-37 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000
IP 0.5-4.0 ug für IP und 1:200-1:1000
für WB
IHC 1:50-1:500

Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : HeLa-Zellen, Mausherzgewebe

IP : HeLa-Zellen,

IHC : humanes Pankreasgewebe,

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser

Hongwei Chu

Pubmed ID

34932315

Journal

Anal Chem

Anwendung

WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HeLa cells were subjected to SDS PAGE followed by western blot with 16578-1-AP (CTSZ antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



IP Result of anti-CTSZ (IP:16578-1-AP, 4ug; Detection:16578-1-AP 1:300) with HeLa cells lysate 2500ug.



Immunohistochemical analysis of paraffin-embedded human pancreas tissue slide using 16578-1-AP (CTSZ antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human pancreas tissue slide using 16578-1-AP (CTSZ antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).