

Nur für Forschungszwecke

HEXIM2 Polyklonaler Antikörper

Katalog-Nr.:16598-1-AP



Allgemeine Informationen

Katalog-Nr.: 16598-1-AP	GenBank-Zugangsnummer: BC025970	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 900 µg/ml von Nanodrop und 400 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 124790	Empfohlene Verdünnungen: WB 1:500-1:1000 IHC 1:20-1:200 IF 1:20-1:200
Wirt: Kaninchen	Vollständiger Name: hexamthylene bis-acetamide inducible 2	
Isotyp: IgG	Berechnete Masse: 286 aa, 32 kDa	
Immunogen Katalognummer: AG9754	Beobachtete Masse: 38 kDa	

Anwendungen

Geprüfte Anwendungen: IF, IHC, WB,ELISA	Positivkontrollen: WB : C6-Zellen, COLO 320-Zellen, HEK-293-Zellen IHC : humanes Hodengewebe, IF : A549-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	
Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



C6 cells were subjected to SDS PAGE followed by western blot with 16598-1-AP (HEXIM2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 16598-1-AP (HEXIM2 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 16598-1-AP (HEXIM2 Antibody) at dilution of 1:50 (under 40x lens).



Immunofluorescent analysis of A549 cells using 16598-1-AP (HEXIM2 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).