

Nur für Forschungszwecke

PRKX Polyklonaler Antikörper

Katalog-Nr.:16651-1-AP



Allgemeine Informationen

Katalog-Nr.:
16651-1-AP

Größe:
150ul, Konzentration: 150 µg/ml von
Nanodrop und 147 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG10083

GenBank-Zugangsnummer:
BC041073

GeneID (NCBI):
5613

Vollständiger Name:
protein kinase, X-linked

Berechnete Masse:
358 aa, 41 kDa

Beobachtete Masse:
41 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

WB: Jurkat-Zellen, Mauslungengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Jurkat cells were subjected to SDS PAGE followed by western blot with 16651-1-AP (PRKX antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.