

Nur für Forschungszwecke

KRT2-specific Polyklonaler Antikörper



Katalog-Nr.:16849-1-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:

16849-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_000423

GeneID (NCBI):

3849

Vollständiger Name:

keratin 2

Berechnete Masse:

65 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus

Zitierte Arten:

Human

Hintergrundinformationen

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are classified into two major sequence types. Type I keratins are a group of acidic intermediate filament proteins, including K9-K23, and the hair keratins Ha1-Ha8. Type II keratins are the basic or neutral counterparts to the acidic type I keratins, including K1-K8, and the hair keratins, Hb1-Hb6. Keratin 2A is a type II cytokeratin. It is found largely in the upper spinous layer of epidermal keratinocytes, and also in other regions, notably the epidermis covering the knee, thigh, and groin. This antibody is specifically against Keratin 2. It does not cross react with other keratins.

Bemerkenswerte Veröffentlichungen

Verfasser

Ningning Dang

Pubmed ID

27096061

Journal

Iran J Basic Med Sci

Anwendung

WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

