

Nur für Forschungszwecke

# KRT12-specific Polyklonaler Antikörper

Katalog-Nr.: 16856-1-AP



## Allgemeine Informationen

Katalog-Nr.:

16856-1-AP

Größe:

150ul, Konzentration: 160 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM\_000223

GeneID (NCBI):

3859

Vollständiger Name:

keratin 12

Berechnete Masse:

54 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

## Hintergrundinformationen

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are classified into two major sequence types. Type I keratins are a group of acidic intermediate filament proteins, including K9–K23, and the hair keratins Ha1–Ha8. Type II keratins are the basic or neutral counterparts to the acidic type I keratins, including K1–K8, and the hair keratins, Hb1–Hb6. Keratin 12 is type I keratin, which is expressed in corneal epithelia. This antibody is specifically against KRT12.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

