

Nur für Forschungszwecke

BTBD10 Polyklonaler Antikörper

Katalog-Nr.:16995-1-AP



Allgemeine Informationen

Katalog-Nr.:
16995-1-AP

Größe:
150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG10684

GenBank-Zugangsnummer:
BC017072

GeneID (NCBI):
84280

Vollständiger Name:
BTB (POZ) domain containing 10

Berechnete Masse:
475 aa, 54 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

BTBD10 is a novel Akt-interacting protein with a BTB/POZ domain, which equally binds to any Akt. BTBD10 appears to behave as a suppressor of cell death including neuronal cell death related to ALS and an enhancer of cell growth via its positive regulation of Akt phosphorylation. The protein is expressed ubiquitously in vivo, distributing in cells as a unique cytoplasmic filamentous around the nucleus. This antibody 16995-1-AP is raised against C-terminal 220 amino acids of full protein.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

