

Nur für Forschungszwecke

# ELP3 Polyklonaler Antikörper

Katalog-Nr.: 17016-1-AP



## Allgemeine Informationen

|                                                            |                                                                        |                                                               |
|------------------------------------------------------------|------------------------------------------------------------------------|---------------------------------------------------------------|
| Katalog-Nr.:<br>17016-1-AP                                 | GenBank-Zugangsnummer:<br>BC001240                                     | Reinigungsmethode:<br>Antigen-Affinitätsreinigung             |
| Größe:<br>150ul, Konzentration: 500 µg/ml von<br>Nanodrop; | GeneID (NCBI):<br>55140                                                | Empfohlene Verdünnungen:<br>WB 1:500-1:2000<br>IHC 1:50-1:500 |
| Wirt:<br>Kaninchen                                         | Vollständiger Name:<br>elongation protein 3 homolog (S.<br>cerevisiae) |                                                               |
| Isotyp:<br>IgG                                             | Berechnete Masse:<br>62 kDa                                            |                                                               |
| Immunogen Katalognummer:<br>AG9993                         | Beobachtete Masse:<br>60-62 kDa                                        |                                                               |

## Anwendungen

|                                                                                                                                                           |                                                                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|
| Geprüfte Anwendungen:<br>IHC, WB, ELISA                                                                                                                   | Positivkontrollen:                                             |
| Getestete Reaktivität:<br>Human, Maus, Ratte                                                                                                              | WB : Maushirngewebe, Rattenhirngewebe<br>IHC : Maushirngewebe, |
| <b>Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.</b> |                                                                |

## Hintergrundinformationen

### Lagerung


Lagerungsbedingungen:  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
Lagerungspuffer:  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 17016-1-AP (ELP3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 17016-1-AP (ELP3 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 17016-1-AP (ELP3 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).