

Nur für Forschungszwecke

APCDD1 Polyklonaler Antikörper

Katalog-Nr.:17024-1-AP



Allgemeine Informationen

Katalog-Nr.:
17024-1-AP

Größe:
150ul, Konzentration: 387 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG10303

GenBank-Zugangsnummer:
BC053324

GeneID (NCBI):
147495

Vollständiger Name:
adenomatosis polyposis coli down-
regulated 1

Berechnete Masse:
514 aa, 59 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
IHC,ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 17024-1-AP (APCDD1 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 17024-1-AP (APCDD1 Antibody) at dilution of 1:50 (under 40x lens).