

Nur für Forschungszwecke

PDLIM7 Polyklonaler Antikörper

Katalog-Nr.:17080-1-AP

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:
17080-1-AP

Größe:
150ul, Konzentration: 300 µg/ml von
Nanodrop und 267 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG1334

GenBank-Zugangsnummer:
BC014521

GeneID (NCBI):
9260

Vollständiger Name:
PDZ and LIM domain 7 (enigma)

Berechnete Masse:
50 kDa

Beobachtete Masse:
48-50 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:2000-1:10000

IP 0.5-4.0 µg für IP und 1:500-1:2000

für WB

IHC 1:20-1:200

Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : K-562-Zellen, HeLa-Zellen

IP : Jurkat-Zellen,

IHC : humanes Lebergewebe, humanes Eierstockgewebe, humanes Hautgewebe, humanes Milzgewebe, humanes Nierengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




WB result of PDLIM7 antibody (17080-1-AP; 1:5000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PDLIM7 transfected HeLa cells.



Various lysates were subjected to SDS PAGE followed by western blot with 17080-1-AP (PDLIM7 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver using 17080-1-AP (PDLIM7 antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver using 17080-1-AP (PDLIM7 antibody) at dilution of 1:100 (under 40x lens).



IP Result of anti-PDLIM7 (IP:17080-1-AP, 4ug; Detection:17080-1-AP 1:1000) with Jurkat cells lysate 3600ug.