

Nur für Forschungszwecke

GBA2 Polyklonaler Antikörper

Katalog-Nr.:17164-1-AP



Allgemeine Informationen

Katalog-Nr.:
17164-1-AP

Größe:
150ul , Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG10606

GenBank-Zugangsnummer:
BC011363

GeneID (NCBI):
57704

Vollständiger Name:
glucosidase, beta (bile acid) 2

Berechnete Masse:
927 aa, 105 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

GBA2(Non-lysosomal glucosylceramidase) is also named as KIAA1605,NLGase and belongs to the non-lysosomal glucosylceramidase family.It catalyzes the conversion of glucosylceramide to free glucose and ceramide and is involved in sphingomyelin generation and prevention of glycolipid accumulation.GBA2 contains a putative membrane-spanning domain and 2 potential N-glycosylation sites. However, glycosidase treatment indicated that GBA2 is not a glycoprotein(PMID:11489889).It has 3 isoforms produced by alternative splicing.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

