

Allgemeine Informationen

Katalog-Nr.: 17213-1-AP	GenBank-Zugangsnummer: BC032008	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 300 µg/ml von Nanodrop;	GeneID (NCBI): 54534	Empfohlene Verdünnungen: WB 1:500-1:2000 IHC 1:50-1:500
Wirt: Kaninchen	Vollständiger Name: mitochondrial ribosomal protein L50	
Isotyp: IgG	Berechnete Masse: 158 aa, 18 kDa	
Immunogen Katalognummer: AG10996	Beobachtete Masse: 18 kDa	

Anwendungen

Geprüfte Anwendungen: IHC, WB, ELISA	Positivkontrollen: WB : HT-1080.Zellen, HeLa-Zellen, HepG2-Zellen IHC : humanes Lebergewebe,
Getestete Reaktivität: Human	

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Hintergrundinformationen


MRPL50, also named as L50mt, is a mitochondrial protein with MW 18 kDa.

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 17213-1-AP (MRPL50 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 17213-1-AP (MRPL50 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).