

Nur für Forschungszwecke

NLN Polyklonaler Antikörper

Katalog-Nr.:17306-1-AP



Allgemeine Informationen

Katalog-Nr.: 17306-1-AP	GenBank-Zugangsnummer: BC004985	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul	GeneID (NCBI): 57486	
Wirt: Kaninchen	Vollständiger Name: neurolysin (metallopeptidase M3 family)	
Isotyp: IgG	Berechnete Masse: 607 aa, 70 kDa	
Immunogen Katalognummer: AG11193		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

NLN(Neurolysin, mitochondrial) is also named as AGTBP, KIAA1226 and belongs to the peptidase M3 family. It is an enzyme involved in the metabolic inactivation of bioactive peptides and functions as a monomer of molecular mass 78 kDa. It also can be detected the molecular weight between 70 kDa and 85 kDa according to the reports(PMID:10095765, 11092934). The full length protein has a transit peptide with 37 amino acids.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

