

UBE2U Polyklonaler Antikörper

Katalog-Nr.: 17323-1-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:
17323-1-AP

Größe:

150ul, Konzentration: 350 µg/ml von
Nanodrop und 373 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG11302

GenBank-Zugangsnummer:

BC029895

GeneID (NCBI):

148581

Vollständiger Name:

ubiquitin-conjugating enzyme E2U
(putative)

Berechnete Masse:

226aa, 27 kDa; 321aa, 38 kDa

Beobachtete Masse:

38 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus

Hintergrundinformationen

UBE2U (Ubiquitin-conjugating enzyme E2 U) is also named as Ubiquitin carrier protein U, Ubiquitin-protein ligase U and belongs to the ubiquitin-conjugating enzyme family. It catalyzes the covalent attachment of ubiquitin to other proteins. It has 2 isoforms produced by alternative splicing with the molecular weight of 38 kDa and 27 kDa.

Bemerkenswerte Veröffentlichungen

Verfasser

Pubmed ID

Journal

Anwendung

Yingying Guo

27903633

J Biol Chem

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

