

Nur für Forschungszwecke

THNSL2 Polyklonaler Antikörper

Katalog-Nr.:17573-1-AP



Allgemeine Informationen

Katalog-Nr.: 17573-1-AP	GenBank-Zugangsnummer: BC064423	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 800 µg/ml von Nanodrop;	GeneID (NCBI): 55258	Empfohlene Verdünnungen: WB 1:500-1:2000
Wirt: Kaninchen	Vollständiger Name: threonine synthase-like 2 (S. cerevisiae)	
Isotyp: IgG	Berechnete Masse: 405 aa, 45 kDa	
Immunogen Katalognummer: AG11788	Beobachtete Masse: 55-60 kDa	

Anwendungen

Geprüfte Anwendungen: WB,ELISA	Positivkontrollen: WB : HeLa-Zellen, HEK-293-Zellen, HepG2-Zellen
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 17573-1-AP (THNSL2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.