

Nur für Forschungszwecke

TSPYL4 Polyklonaler Antikörper

Katalog-Nr.:17576-1-AP



Allgemeine Informationen

Katalog-Nr.: 17576-1-AP	GenBank-Zugangsnummer: BC009116	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 650 µg/ml von Nanodrop;	GeneID (NCBI): 23270	Empfohlene Verdünnungen: WB 1:2000-1:10000
Wirt: Kaninchen	Vollständiger Name: TSPY-like 4	
Isotyp: IgG	Berechnete Masse: 414 aa, 45 kDa	
Immunogen Katalognummer: AG10730	Beobachtete Masse: 45 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Maus, Ratte	WB : Maus-Cerebellum-Gewebe, Ratten-Cerebellum-Gewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 17576-1-AP (TSPYL4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.