

Nur für Forschungszwecke

GRHL1 Polyklonaler Antikörper

Katalog-Nr.:17644-1-AP



Allgemeine Informationen

Katalog-Nr.:
17644-1-AP

Größe:
150ul, Konzentration: 200 µg/ml von
Nanodrop und 167 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG11832

GenBank-Zugangsnummer:
BC067519

GeneID (NCBI):
29841

Vollständiger Name:
grainyhead-like 1 (Drosophila)

Berechnete Masse:
618 aa, 70 kDa

Beobachtete Masse:
70 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

WB: NCCIT cell,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



NCCIT cell were subjected to SDS PAGE followed by western blot with 17644-1-AP (GRHL1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.