

Nur für Forschungszwecke

# Jagged1 Polyklonaler Antikörper

Katalog-Nr.:17791-1-AP



## Allgemeine Informationen

<b>Katalog-Nr.:</b> 17791-1-AP	<b>GenBank-Zugangsnummer:</b> BC098393	<b>Reinigungsmethode:</b> Antigen-Affinitätsreinigung
<b>Größe:</b> 150ul , Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	<b>GeneID (NCBI):</b> 182	
<b>Wirt:</b> Kaninchen	<b>Vollständiger Name:</b> jagged 1 (Alagille syndrome)	
<b>Isotyp:</b> IgG	<b>Berechnete Masse:</b> 1218 aa, 134 kDa	
<b>Immunogen Katalognummer:</b> AG11891		

## Anwendungen

**Geprüfte Anwendungen:**  
ELISA

**Getestete Reaktivität:**  
Human, Maus, Ratte

## Hintergrundinformationen

JAG1, also named as JAGL1, CD339 and Jagged1, is a ligand for multiple Notch receptors and involved in the mediation of Notch signaling. It may be involved in cell-fate decisions during hematopoiesis. JAG1 seems to be involved in early and late stages of mammalian cardiovascular development. It inhibits myoblast differentiation. JAG1 enhances fibroblast growth factor-induced angiogenesis. JAG1 signalling through notch 1 has also been shown to play a role in hematopoiesis.

## Lagerung

**Lagerungsbedingungen:**  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

**Lagerungspuffer:**  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

