

Nur für Forschungszwecke

KCNN1 Polyklonaler Antikörper

Katalog-Nr.:17929-1-AP



Allgemeine Informationen

Katalog-Nr.:
17929-1-AP

Größe:
150ul, Konzentration: 700 µg/ml von
Nanodrop;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG12357

GenBank-Zugangsnummer:
BC075037

GeneID (NCBI):
3780

Vollständiger Name:
potassium intermediate/small
conductance calcium-activated
channel, subfamily N, member 1

Berechnete Masse:
543 aa, 60 kDa

Beobachtete Masse:
70 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

WB: Maushirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse brain tissue were subjected to SDS PAGE followed by western blot with 17929-1-AP (KCNN1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.