

Nur für Forschungszwecke

TNFAIP1 Polyklonaler Antikörper

Katalog-Nr.: 18234-1-AP



Allgemeine Informationen

Katalog-Nr.:

18234-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG12941

GenBank-Zugangsnummer:

BC003694

GeneID (NCBI):

7126

Vollständiger Name:

tumor necrosis factor, alpha-induced protein 1 (endothelial)

Berechnete Masse:

36 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

TNFAIP1 was identified as a human gene activated transcriptionally by tumour necrosis factor alpha in endothelial cells. It is encoded by a 3.5-kilobase transcript and its expression is found to be developmentally regulated in a tissue-specific manner. The open reading frame predicts a 316-amino acid sequence rich in charged residues, particularly at the carboxyl terminus. The expected band is about 36-40kd.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

