

Nur für Forschungszwecke

SPEG Polyklonaler Antikörper

Katalog-Nr.:18716-1-AP



Allgemeine Informationen

Katalog-Nr.:
18716-1-AP

Größe:
150ul

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_005876

GeneID (NCBI):
10290

Vollständiger Name:
SPEG complex locus

Berechnete Masse:
354 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

SPEG, also named as BPEG, APEG1 and KIAA1297, is highly similar to the corresponding rat and mouse proteins. Expression of SPEG is thought to serve as a marker for differentiated vascular smooth muscle cells which may have a role in regulating growth and differentiation of this cell type. SPEG is expressed in vascular smooth muscle cells and is down-regulated by vascular injury. Multiple transcripts are produced from the SPEG complex locus by alternative splicing and alternative promoter usage. The antibody recognizes the C-term of SPEG.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

