

Nur für Forschungszwecke

HOXD11 Polyklonaler Antikörper

Katalog-Nr.:18734-1-AP

2 Publikationen



Allgemeine Informationen

Katalog-Nr.:
18734-1-AP

Größe:
150ul, Konzentration: 450 µg/ml von
Nanodrop;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
BC160135

GeneID (NCBI):
3237

Vollständiger Name:
homeobox D11

Berechnete Masse:
35 kDa

Beobachtete Masse:
32-36 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

In Publikationen genannte Anwendungen:
IF, IHC

Getestete Reaktivität:
Human, Maus, Ratte

Zitierte Arten:
Human

Positivkontrollen:

WB : SMMC-7721-Zellen,

Hintergrundinformationen

HOXD11, also named as HOX4F, belongs to the Abd-B homeobox family. HOXD11 is a sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. HOXD11 plays a role in forelimb morphogenesis.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Karim Fekir	30837223	Cancer Res	IHC
Jian Hui Low	31303547	Cell Stem Cell	IF

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 18734-1-AP (HOXD11 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.