

Nur für Forschungszwecke

PROP1 Polyklonaler Antikörper

Katalog-Nr.:18739-1-AP



Allgemeine Informationen

Katalog-Nr.:

18739-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

BC069076

GeneID (NCBI):

5626

Vollständiger Name:

PROP paired-like homeobox 1

Berechnete Masse:

25 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

PROP1 belongs to the paired homeobox family. It contains one homeobox DNA-binding domain. PROP1 is possibly involved in the ontogenesis of pituitary gonadotropes, as well as somatotropes, lactotropes and caudomedial thyrotropes. It is expressed specifically in embryonic pituitary. Mutation of PROP1 will cause combined pituitary hormone deficiency (CPHD) which is characterized by impaired production of growth hormone (GH) and one or more of the other five anterior pituitary hormones.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

