

Nur für Forschungszwecke

EGFL7 Polyklonaler Antikörper

Katalog-Nr.: 19291-1-AP

Vorgestelltes Produkt

11 Publikationen



Allgemeine Informationen

Katalog-Nr.:

19291-1-AP

Größe:

150ul, Konzentration: 167 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG6559

GenBank-Zugangsnummer:

BC012377

GeneID (NCBI):

51162

Vollständiger Name:

EGF-like-domain, multiple 7

Berechnete Masse:

273 aa, 30 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

IHC, ELISA

In Publikationen genannte Anwendungen:

IF, IHC, WB

Getestete Reaktivität:

Human

Zitierte Arten:

Human, Maus, Ratte

Hintergrundinformationen

Epidermal growth factor-like domain-containing protein 7 (EGFL7), also known as vascular endothelial statin, is an endothelial cell-derived secreted factor that regulates vascular tube formation. This protein may be involved in the growth and proliferation of tumor cells. Recent studies have reported elevated expression of EGFL7 in several tumors and cancer cell lines, including kidney tumors, malignant gliomas, hepatocellular carcinomas, and colon cancers.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Xiao-Hong Wang	25270395	Int J Mol Med	WB
Caterina L Mammola	30504933	Eur J Histochem	IHC
Chunqing Yang	32464549	Mol Ther Nucleic Acids	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 19291-1-AP (EGFL7 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 19291-1-AP (EGFL7 Antibody) at dilution of 1:50 (under 40x lens).