

Nur für Forschungszwecke

TLR1-Specific Polyklonaler Antikörper

Katalog-Nr.: 19500-1-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:

19500-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_003263

GeneID (NCBI):

7096

Vollständiger Name:

toll-like receptor 1

Berechnete Masse:

90 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Zitierte Arten:

Ratte

Hintergrundinformationen

TLR1, also named as TIL, CD281 and KIAA0012, belongs to the Toll-like receptor family. TLR1 participates in the innate immune response to microbial agents. TLR1 cooperates with TLR2 mediates the innate immune response to bacterial lipoproteins or lipopeptides. Acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response. This antibody is specific to TLR1.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Dongmin Li	25534427	Mol Cell Endocrinol	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

