

Nur für Forschungszwecke

# PTPN20B Polyklonaler Antikörper

Katalog-Nr.:19510-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
19510-1-AP

Größe:  
150ul

Wirt:  
Kaninchen

Isotyp:  
IgG

GenBank-Zugangsnummer:  
NM\_015605

GeneID (NCBI):  
26095

Vollständiger Name:  
protein tyrosine phosphatase, non-receptor type 20B

Berechnete Masse:  
48 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human

## Hintergrundinformationen

PTPN20B, also named as PTPN20A, belongs to the protein-tyrosine phosphatase family and Non-receptor class subfamily. PTPN20B is a tyrosine-protein phosphatase targeted to sites of actin polymerization in response of varied extracellular stimuli. PTPN20B has tyrosine phosphatase activity towards various tyrosyl phosphorylated substrates. PTPN20B catalyze the reaction: Protein tyrosine phosphate + H<sub>2</sub>O = protein tyrosine + phosphate. The antibody can recognize the isoform 1,2,7,8,10,11,13,15.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

