

Nur für Forschungszwecke

# VPS15-phospho(S1079) Polyklonaler Antikörper



Katalog-Nr.:19552-1-AP

## Allgemeine Informationen

**Katalog-Nr.:**  
19552-1-AP

**Größe:**  
100ul , Konzentration: 493 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

**Wirt:**  
Kaninchen

**Isotyp:**  
IgG

**GenBank-Zugangsnummer:**  
NM\_014602

**GeneID (NCBI):**  
30849

**Vollständiger Name:**  
phosphoinositide-3-kinase,  
regulatory subunit 4

**Berechnete Masse:**  
153 kDa

**Reinigungsmethode:**  
Antigen-Affinitätsreinigung

## Anwendungen

**Geprüfte Anwendungen:**  
ELISA

**Getestete Reaktivität:**  
Human, Maus, Ratte

## Hintergrundinformationen

PIK3R4, also named as PI3-kinase p150 subunit, belongs to the protein kinase superfamily and Ser/Thr protein kinase family. PIK3R4 is regulatory subunit of the PI3K complex. PIK3R4 may regulate membrane trafficking late in the endocytic pathway. PIK3R4 catalyze the reaction: ATP + a protein = ADP + a phosphoprotein. The antibody is specific to Ser1079(phospho).

## Lagerung

**Lagerungsbedingungen:**

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

**Lagerungspuffer:**

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

