

Nur für Forschungszwecke

SDPR-Specific Polyklonaler Antikörper

Katalog-Nr.: 19558-1-AP



Allgemeine Informationen

Katalog-Nr.:

19558-1-AP

Größe:

150ul, Konzentration: 253 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_004657

GeneID (NCBI):

8436

Vollständiger Name:

serum deprivation response
(phosphatidylserine binding protein)

Berechnete Masse:

47 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Serum deprivation-response protein (SDRP), a substrate of protein kinase C, is a member of a group of proteins that bind to phosphatidylserine in a calcium-dependent manner. SDRP may participate in the generation or stabilization of membrane curvature. SDRP is highly expressed in heart and lung, and low in brain, kidney, liver, pancreas, placenta and skeletal muscle. This is a rabbit polyclonal anti-SDRP antibody raised against C-terminus of human SDRP.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

