

Nur für Forschungszwecke

INCENP Polyklonaler Antikörper

Katalog-Nr.:19689-1-AP



Allgemeine Informationen

Katalog-Nr.: 19689-1-AP	GenBank-Zugangsnummer: NM_001040694	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul	GeneID (NCBI): 3619	
Wirt: Kaninchen	Vollständiger Name: inner centromere protein antigens	
Isotyp: IgG	Berechnete Masse: 105 kDa	

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

INCENP belongs to the INCENP family. It's component of the chromosomal passenger complex (CPC), a complex that acts as a key regulator of mitosis. The CPC complex has essential functions at the centromere in ensuring correct chromosome alignment and segregation and is required for chromatin-induced microtubule stabilization and spindle assembly. INCENP probably acts through association with AURKB or AURKC. It seems to bind directly to microtubules. The antibody is to INCENP with Y(822-phospho) and S(831-phospho).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

