

Allgemeine Informationen

Katalog-Nr.:

19693-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_017629

GeneID (NCBI):

192670

Vollständiger Name:

eukaryotic translation initiation

factor 2C, 4

Berechnete Masse:

97 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

Hintergrundinformationen

EIF2C4, also named as AGO4 and KIAA1567, belongs to the argonaute family and Ago subfamily. It is required for RNA-mediated gene silencing (RNAi). EIF2C4 binds to short RNAs such as microRNAs (miRNAs) and represses the translation of mRNAs which are complementary to them. EIF2C4 lacks endonuclease activity and does not appear to cleave target mRNAs. It is also required for RNA-directed transcription and replication of the human hepatitis delta virus (HDV).

Bemerkenswerte Veröffentlichungen

Verfasser

Pubmed ID

Journal

Anwendung

Lianjiu Su

31694915

J Biol Chem

WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

