

Nur für Forschungszwecke

NFIX-Specific Polyklonaler Antikörper



Katalog-Nr.: 19708-1-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:

19708-1-AP

Größe:

150ul, Konzentration: 933 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_002501

GeneID (NCBI):

4784

Vollständiger Name:

nuclear factor I/X (CCAAT-binding
transcription factor)

Berechnete Masse:

55 kDa

Beobachtete Masse:

23 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

Hintergrundinformationen

NFIX, also named as CTF, belongs to the CTF/NF-I family. NFIX recognizes and binds the palindromic sequence 5'-TTGGCCNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. NFIX is individually capable of activating transcription and replication. The antibody is specific to NFIX.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
I Wayan Sumardika	30720226	Mol Carcinog	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

