

Nur für Forschungszwecke

PPP1R3A-Specific Polyklonaler Antikörper

Katalog-Nr.: 19712-1-AP



Allgemeine Informationen

Katalog-Nr.:

19712-1-AP

Größe:

150ul, Konzentration: 167 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_002711

GeneID (NCBI):

5506

Vollständiger Name:

protein phosphatase 1, regulatory
(inhibitor) subunit 3A

Berechnete Masse:

126 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

PPP1R3A, also named as PP1G, seems to act as a glycogen-targeting subunit for PP1. PP1 is essential for cell division, and participates in the regulation of glycogen metabolism, muscle contractility and protein synthesis. PPP1R3A plays an important role in glycogen synthesis but is not essential for ins activation of glycogen synthase. This antibody is specific to PP1R3A.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

