

Nur für Forschungszwecke

RASGRF2-Specific Polyklonaler Antikörper



Katalog-Nr.: 19717-1-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:

19717-1-AP

Größe:

150ul, Konzentration: 140 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_006909

GeneID (NCBI):

5924

Vollständiger Name:

Ras protein-specific guanine
nucleotide-releasing factor 2

Berechnete Masse:

141 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

RASGRF2, also named as GRF2, functions as a calcium-regulated nucleotide exchange factor activating both Ras and RAC1 through the exchange of bound GDP for GTP. RASGRF2 preferentially activates HRAS in vivo compared to RRAS based on their different types of prenylation. RASGRF2 functions in synaptic plasticity by contributing to the induction of long term potentiation. This antibody is specific to RASGRF2.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Kari A Herrington	28539409	Mol Biol Cell	

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

