

Nur für Forschungszwecke

# SMARCA2-Specific Polyklonaler Antikörper



Katalog-Nr.:19718-1-AP

## Allgemeine Informationen

Katalog-Nr.:

19718-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM\_003070

GeneID (NCBI):

6595

Vollständiger Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

Berechnete Masse:

181 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

## Hintergrundinformationen

SMARCA2 is a transcriptional coactivator cooperating with nuclear hormone receptors to potentiate transcriptional activation. SMARCA2 is involved in vitamin D-coupled transcription regulation via its association with the WINAC complex, a chromatin-remodeling complex recruited by vitamin D receptor (VDR), which is required for the ligand-bound VDR-mediated transrepression of the CYP27B1 gene. The antibody is specific to SMARCA2

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

