

Nur für Forschungszwecke

SNTB1-Specific Polyklonaler Antikörper

Katalog-Nr.:19719-1-AP



Allgemeine Informationen

Katalog-Nr.:
19719-1-AP

Größe:
150ul, Konzentration: 287 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_021021

GeneID (NCBI):
6641

Vollständiger Name:
syntrophin, beta 1 (dystrophin-
associated protein A1, 59kDa, basic
component 1)

Berechnete Masse:
58 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

SNTB1, also named as SNT2B1, BSYN2 and TIP-43, belongs to the syntrophin family. SNTB1 is an adapter protein that binds to and probably organizes the subcellular localization of a variety of membrane proteins. It may link various receptors to the actin cytoskeleton and the dystrophin glycoprotein complex. This antibody is specific to SNTB1.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

