

Nur für Forschungszwecke

TAF4-Specific Polyklonaler Antikörper

Katalog-Nr.:19725-1-AP



Allgemeine Informationen

Katalog-Nr.:
19725-1-AP

Größe:
150ul, Konzentration: 140 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_003185

GeneID (NCBI):
6874

Vollständiger Name:
TAF4 RNA polymerase II, TATA box
binding protein (TBP)-associated
factor, 135kDa

Berechnete Masse:
110 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

TAF4, also named as TAF2C, TAF2C1, TAF4A, TAFII130 and TAFII135, belongs to the TAF4 family. TAF4 makes part of TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. TAF4 potentiates transcriptional activation by the AF-2S of the retinoic acid, vitamin D3 and thyroid hormone. This antibody is specific to TAF4.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

