

Nur für Forschungszwecke

TAF5-Specific Polyklonaler Antikörper

Katalog-Nr.:19726-1-AP



Allgemeine Informationen

Katalog-Nr.: 19726-1-AP	GenBank-Zugangsnummer: NM_006951	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 213 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 6877	
Wirt: Kaninchen	Vollständiger Name: TAF5 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 100kDa	
Isotyp: IgG	Berechnete Masse: 87 kDa	

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

TAF5, also named as TAF2D and TAFII100, belongs to the WD repeat TAF5 family. TAFs are components of the transcription factor IID (TFIID) complex, PCAF histone acetylase complex and TBP-free TAFII complex (TFTC). TAFs components-TFIIFD is essential for mediating regulation of RNA polymerase transcription. TAF5 interacts strongly with the histone H4-related TAF6/TAFII80 and the histone H3-related TAF9/TAFII31, as well as a stable complex comprised of both TAF5/TAFII80 and TAF6/TAFII31. Apparently weaker interactions of TAF5/TAFII100 with TBP, TAF1/TAFII250, TAF11/TAFII28, and TAF12/TAFII20, but not TAF7/TAFII55, also have been observed. This antibody is specific to TAF5.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

