

Nur für Forschungszwecke

# RASGRP2 Polyklonaler Antikörper

Katalog-Nr.:19745-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
19745-1-AP

Größe:  
150ul, Konzentration: 133 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

GenBank-Zugangsnummer:  
NM\_153819

GeneID (NCBI):  
10235

Vollständiger Name:  
RAS guanyl releasing protein 2  
(calcium and DAG-regulated)

Berechnete Masse:  
69 kDa

Beobachtete Masse:  
48 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human, Maus

## Hintergrundinformationen

RASGRP2, also named as CDC25L and MCG7, belongs to the RASGRP family. It functions as a calcium- and DAG-regulated nucleotide exchange factor specifically activating Rap through the exchange of bound GDP for GTP. RASGRP2 may also activate other GTPases such as RRAS, RRAS2, NRAS, KRAS but not HRAS. RASGRP2 functions in aggregation of platelets and adhesion of T-lymphocytes and neutrophils probably through inside-out integrin activation. It may function in the muscarinic acetylcholine receptor M1/CHRM1 signaling pathway. The antibody recognizes the C-term of RASGRP2.

## Lagerung

Lagerungsbedingungen:  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
Lagerungspuffer:  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

