

Nur für Forschungszwecke

Collagen Type XIII Polyklonaler Antikörper



Katalog-Nr.: 19801-1-AP

Allgemeine Informationen

Katalog-Nr.:

19801-1-AP

Größe:

150ul, Konzentration: 467 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_001130103

GeneID (NCBI):

1305

Vollständiger Name:

collagen, type XIII, alpha 1

Berechnete Masse:

70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

COL13A1 is involved in cell-matrix and cell-cell adhesion interactions that are required for normal development. COL13A1 may participate in the linkage between muscle fiber and basement membrane. COL13A1 may play a role in endochondral ossification of bone and branching morphogenesis of lung. COL13A1 can bind heparin. The antibody can recognize all the isoforms except isoform9 of COL13A1. This antibody is specific to COL13A1.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

